

How to run the Online Diagnostic Lab Management System Project (ODLMS)

1. Open the zip file
2. Extract the file and copy odlms folder(inside the source_code)
3. Paste inside root directory(for xampp xampp/htdocs, for wamp wamp/www, for lamp var/www/html)
4. Open PHPMyAdmin (<http://localhost/phpmyadmin>)
5. Create a database with name odlmsdb
6. Import odlmsdb.sql file(given inside the zip package in SQL file folder)
7. Run the script <http://localhost/odlms>

Admin Credential

Username: dbms_mini_project

Password:dbms12345

Employee Credential

Username: lab001

Password: abhi123

User/ pateint Credential

Username: drupad1522@gmail.com

Password: drupad1522